


Plan Quality Algorithm: RTOG 0529 N+ Anal Canal [35 Max Possible] [35 Metrics] (Page 1 of 4)












Plan Quality Algorithm: RTOG 0529 N+ Anal Canal [35 Max Possible] [35 Metrics] (Page 2 of 4)

11	 V (%) D(Gy)	Volume (%) of the PTV45 covered by 45 (Gy)	1	
12	 V (%) D(Gy)	Volume (%) of the PTV45 covered by 40.5 (Gy)	1	
13	 V (%) D(Gy)	Volume (%) of the PTV45 covered by 36 (Gy)	1	
14	 V (cc) D(Gy)	Volume (cc) of the SMALL BOWEL covered by 30 (Gy)	1	
15	 V (cc) D(Gy)	Volume (cc) of the SMALL BOWEL covered by 35 (Gy)	1	
16	 V (cc) D(Gy)	Volume (cc) of the SMALL BOWEL covered by 45 (Gy)	1	
17	 ROI Max (Gy)	Maximum dose (Gy) inside the SMALL BOWEL	1	
18	 V (%) D(Gy)	Volume (%) of the RIGHT FEMORAL HEAD covered by 30 (Gy)	1	
19	 V (%) D(Gy)	Volume (%) of the RIGHT FEMORAL HEAD covered by 40 (Gy)	1	
20	 V (%) D(Gy)	Volume (%) of the RIGHT FEMORAL HEAD covered by 44 (Gy)	1	

Plan Quality Algorithm: RTOG 0529 N+ Anal Canal [35 Max Possible] [35 Metrics] (Page 3 of 4)

21	V (%) D(Gy)	Volume (%) of the LEFT FEMORAL HEAD covered by 30 (Gy)	1	≤ 50 ≥ 50.00
22	V (%) D(Gy)	Volume (%) of the LEFT FEMORAL HEAD covered by 40 (Gy)	1	≤ 35 ≥ 35.00
23	V (%) D(Gy)	Volume (%) of the LEFT FEMORAL HEAD covered by 44 (Gy)	1	≤ 5 < 10 ≥ 10
24	V (%) D(Gy)	Volume (%) of the ILIAC CRESTS covered by 30 (Gy)	1	≤ 50 ≥ 50.00
25	V (%) D(Gy)	Volume (%) of the ILIAC CRESTS covered by 40 (Gy)	1	≤ 35 ≥ 35.00
26	V (%) D(Gy)	Volume (%) of the ILIAC CRESTS covered by 50 (Gy)	1	≤ 5 ≥ 5.001
27	V (%) D(Gy)	Volume (%) of the EXTERNAL GENITALIA covered by 20 (Gy)	1	≤ 50 ≥ 50.00
28	V (%) D(Gy)	Volume (%) of the EXTERNAL GENITALIA covered by 30 (Gy)	1	≤ 35 ≥ 35.00
29	V (%) D(Gy)	Volume (%) of the EXTERNAL GENITALIA covered by 40 (Gy)	1	≤ 5 ≥ 5.001
30	V (%) D(Gy)	Volume (%) of the BLADDER covered by 35 (Gy)	1	≤ 50 ≥ 50.00

Plan Quality Algorithm: RTOG 0529 N+ Anal Canal [35 Max Possible] [35 Metrics] (Page 4 of 4)

31	 V (Gy) D (Gy)	Volume (%) of the BLADDER covered by 40 (Gy)	1	
32	 V (Gy) D (Gy)	Volume (%) of the BLADDER covered by 50 (Gy)	1	
33	 V (cc) D (Gy)	Volume (cc) of the LARGE BOWEL covered by 30 (Gy)	1	
34	 V (cc) D (Gy)	Volume (cc) of the LARGE BOWEL covered by 35 (Gy)	1	
35	 V (cc) D (Gy)	Volume (cc) of the LARGE BOWEL covered by 45 (Gy)	1	